For more information
on this product please
scan QR code

Instructions For Use

**Návod k použití • Gebrauchsanweisung • Instrucciones de uso
Mode d'emploi • Οδηγίες χρήσης • Istruzioni per l'uso
Gebruiksaanwijzing • Instruksjoner for bruk • Instrukcja użycia
Instruções de uso • Användningsinstruktioner**

Disposable Veress Needle

DE	Einweg-Veress-Nadel
ES	Aguja de Veress desechable
FR	Aiguille de Veress jetable
IT	Ago di Veress monouso
NL	Wegwerpbare Veress-naald
NO	Engangs Veress-nål
PT	Agulha de Veress descartável
SE	Veress-nål för engångsbruk

Indications / Intended use

Fairmont's Disposable Veress Needle is intended for use in minimally invasive surgical procedures, predominantly abdominal laparoscopic techniques. It is used to introduce or remove gas from the peritoneal cavity as a therapeutic or surgical/radiological procedural method. It is often inserted into the peritoneal cavity for the purpose of insufflation (e.g., with CO₂) to establish pneumoperitoneum prior to abdominal endoscopy.

Instructions For Use

1. Peel open pouch and remove product using standard aseptic techniques.
2. Prior to insertion, inspect the Veress Needle to ensure that the cannula moves freely over of the blunt obturator with minimal pressure applied to the device.
3. Insert the Veress Needle.
4. It is recommended to perform various checks to ensure that the Veress Needle is in the abdominal cavity ensuring red indicator is **not** visible on the safety window.
5. Once correct placement is ascertained, inflate to desired pressure.
6. Once insufflation is complete and the desired abdominal pressure is accomplished, remove the Veress Needle and proceed with laparoscopic procedure.


Precautions

- Single use only, do not re-use, re-process or re-sterilise.
- Suitably trained healthcare personnel may only use this device.
- Do not use if the stylet does not return over the blunt end when pressure on the stylet is released.
- Do not use the Veress Needle if the stylet is not functioning properly.
- This device is sterile (ETO). Please note, this IFU pertains to sterile devices sold by Fairmont Medical. If the device has code "DVNxxxNS" please refer to 3rd party manufacturer for sterility specifications.
- This device must be disposed of in accordance with local environmental laws and regulations

Materials

PVC, Acrylic, Natural TPR, stainless steel.

All components and packaging are latex and Phthalate-free PVC.

Any serious incidents that occur in relation to the device should be reported to the manufacturer and the competent authority of the Member State of the patient and/or user.

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
Is a medical device

EU Importer

Authorised representative in the European Community

Product conforms with Regulation (EU) 2017/745

Language denotation as per EC language codes

Do not use if packaging is damaged

Caution is necessary

Keep away from sunlight

Keep dry

Single use only

Do not resterilise

Non-sterile

Consult instructions for use


Made in Australia

Latex free

Phthalate and DEHP free

Upper and lower temperature limits

Non-pyrogenic

Single sterile barrier system with protective packaging inside, sterilised with EO

 Fairmont Medical Products Pty Ltd
 90 Prosperity Way
 Dandenong South
 Victoria 3175
 Australia
 Ph +61 3 9720 8840
 www.fairmontmedical.com.au

Fairmont Medical Products Ltd
 Elizabeth House
 28 Baddow Road Chelmsford
 Ph +44 (0) 1245 206800
 www.fairmontmedical.com

Fairmont Medical Products Pty Ltd is certified to ISO13485:2016

PM069 Version 4.0

Issued by L. Wood 16/11/2023

Product specifications subject to change without notice. Pictures shown may vary from actual product.
© 2018 Fairmont Medical Products Pty Ltd all rights reserved.