




PRODUCT	FAIRMONT CODE	NPC CODE	PRODUCT DESCRIPTION	BOX QTY	SUITABILITY	SIZE
---------	---------------	----------	---------------------	---------	-------------	------

CAMERA COVERS

	DCC8007	FBE044	Camera Cover Sterile 100% leak free enclosed adapter with Telescopic Lock	40	All Operating Camera Systems	2850mm x140mm
	DCC8000	VEB336	Perforated Tip(Telescopically Folded)	40	All Operating Camera Systems	1850mmx140mm
	DCC8005	VEB337	Enclosed Tip(Telescopically Folded) with vinyl end and Lens	40	All Operating Camera Systems	2100mmx140mm and vinyl sock 220mm
	DCC8006	VEB338	Enclosed Tip(Telescopically Folded) with vinyl end without Lens	40	All Operating Camera Systems	2100mmx140mm and vinyl sock 220mm







DENTAL DRILL COVERS

	DDC8001	VEB339	Perforated Telescopically Folded	40		1800mm length 60mm width
	DDC8003	VEB340	Non Perforated Telescopically Folded	40		1800mm length 60mm width


X-RAY COVERS

	DXC8015	VEB341	Universal	20		1500mm length 1040mm width
	DXC8025	VEB342	Shower Cap style	50		500mm length 1040mm width
	DXC8027	VEB343	Shower Cap style with securing ribbon	20		800mm length 1040mm width





CO2 INSUFFLATOR FILTERS & TUBING SETS

	DIF7002	FTK6240	30mm Silicone Universal Adapter	50	6mm Barb	Tubing Length: 3m with rotating Male Luer Lock
	DIF7004	FTK6236	22mm Female ISO	50	22mm Male Outlet	Tubing Length: 3m with rotating Male Luer Lock
	DIF7005	FTK6244	Male Luer Lock Connector	50	Female Luer Outlet	Tubing Length: 3m with rotating Male Luer Lock
	DIF7007	FTK6246	Male Luer Lock Connector Connector with 1 white hexagonal connector	50	Sopro Outlet	Tubing Length: Y Extension & 3m Tubing with rotating Male Luer Lock
	DIF7009	FTK042	Female Luer Lock	45	Male Luer Lock	Tubing Length: 3m with rotating Male Luer Lock
	DIF7205	FTK6235	15mm Male Connector	50	15mm Female Outlets	Tubing Length: 3m with rotating Male Luer Lock

VERESS NEEDLE

	DVN7012	FTK6267	Colour Grip: Blue	30		Length: 120mm
	DVN7015	FTK6268	Colour Grip: Yellow	30		Length: 150mm
	DVN7200	FTK6266	Colour Grip: Orange	30		Length: 200mm

LAPAROSCOPIC IRRIGATION SUCTION SETS

	LISS4004	FTK6238	Single Use Irrigation Suction Probe with White Connector	5		Length: 380mm Probe OD: 3mm
	LISS4017	FTK6239	Single Use Irrigation Suction Probe with White Connector	10		Length: 440mm Probe OD: 5mm
	LISS4117	FTK6237	Single Use Irrigation Suction Probe with White Connector	5		Length: 380mm Probe OD: 10mm
	LISS4121	FTK6243	LISS with detached 5mm probe	10		Length: 3.7m Single Spike
	LISS4122	FTK6242	LISS with detached 5mm probe	10		Length: 3.7m Double Spike
	LISS4422	FTK6241	LISS with detached 10mm probe	10		Length: 3.7m Double Spike



Fairmont
MEDICAL

fairmontmedical.com.au

Corporate Office
Fairmont Medical
11 Scoresby Road
Bayswater 3153 VIC
Ph 1300 972 088
Ph International +613 9720 8840
Fx International +613 9720 8860

UK Head Office
Fairmont Medical Products Ltd
Elizabeth House
28 Baddow Road Chelmsford
Essex CM2 0DG
Ph: +44 (0)1245 206800

