

# Diathermy Smoke Pen

Fairmont Medical's Diathermy Smoke Pen optimises smoke evacuation contributing to a healthier O.R.

Fairmont Medical offers an economic alternative smoke evacuation pen. **Fixie** provides the opportunity to utilise a safe and efficient method of surgical diathermy smoke removal. Available in 3 size nozzle lengths, **Fixie** provides choice for those seeking a cost effective solution.

- Choice of 3 interchangeable transparent smoke evacuation nozzles
- A selection of loop and ball hooded electrodes available
- Super wide smoke evacuation coverage > 50mm
- Aerodynamic design affords high vacuum > 80L/min
- Light, flexible tubing and swivel connector enables enhanced handling and control
- Easy to identify colour coded cut and coag buttons
- Compact streamlined fully integrated hose and cable
- Saving valuable operating time
- Single use, sterile, ready to use
- Phthalate (DEHP) and latex free
- Australian designed quality



**Fixie Pen**

**Interchangeable Nozzles**

**Hooded Electrodes**



Pencavac Diathermy Smoke FIXIE Pen - Sterile, Disposable			
Product	Product Code	Description	Box Qty
	DSP2130	Disposable Smoke Evacuator Diathermy Pen with 70mm Teflon Blade electrode, Smoke Evacuation Nozzle, 3.0m Tubing with 22mm Female/Male connector and 3m cable with 3-prong plug	10
	DSP2830	Disposable Smoke Evacuator Diathermy Pen with 70mm + 120mm + 152mm Teflon Blades electrode, Smoke Evacuation Nozzles, 3.0m Tubing with 22mm Female/Male connector and 3m cable with 3-prong plug	10